

DISPOSABLES



HUMAN SEMEN ANALYSIS

SCASCOPE

www.scascope.com



MICROPTIC
AUTOMATIC DIAGNOSTIC SYSTEMS

mot

MOTILITY AND CONCENTRATION

Counting chambers SCA

Disposable counting chambers to analyse motility, concentration and pH of semen:

- 10 micron – 2 chambers – pH
- 10 micron – 4 chambers
- 10 micron – 6 chambers
- 20 micron – 2 chambers
- 20 micron – 6 chambers



References:	MOT-10-2-PH
MOT-10-4	MOT-10-6
MOT-20-2	MOT-20-6

mrf

MORPHOLOGY

Morphology slides SCA

Pre-stained slides to analyse morphology of sperm that combine conventional microscope slides and isosmotic SpermBlue stain into a ready to use device.



Reference:	MRF-SB
------------	--------

Use of mounting medium required (Reference: EUK-NEO-100).

vit

VITALITY

Vitality slides SCA

Eosine-Nigrosine based stain, *BrightVit*, optimised in one-step staining technique for the assessment of sperm vitality in brightfield.



Reference:	VIT-BV
------------	--------

dna

DNA FRAGMENTATION

GoldCyto DNA

Test to assess sperm DNA fragmentation, based on the Sperm Chromatin Dispersion (SCD) test. 20 analysis.



Reference:	GC-SCD-20
------------	-----------

Not available in some countries due to patent reasons.

leu

LEUKOCYTES

Leukocytes slides SCA

Kit for detection of peroxidase positive cells in semen. A semen sample is mixed with a substrate specific for the enzyme peroxidase.



Reference:	LEU-SLIDES
------------	------------